

# Technical Data Sheet

## InVivoMAb anti-human/rat/fish AChR



**Attention:** Use of this product constitutes an agreement to Bio X Cell's Terms and Conditions which are included with this product in print and can also be found at <https://bioxcell.com/terms-and-conditions>.

### Lot Specific Information

Lot Number: Lot Specific\*  
Volume: Lot Specific\*  
Concentration: Lot Specific\* (generally 4 to 11 mg/ml) \*  
Total Protein: Lot Specific\*

\*This information will be noted on the certificate of analysis that ships with this product.

### Product Information

Catalog Number: BE0123  
Clone: Mab35 (TIB-175)  
Isotype: Rat IgG1  
Recommended Isotype Control(s): InVivoMAb rat IgG1 isotype control, anti-horseradish peroxidase  
Recommended Dilution Buffer: InVivoPure pH 7.0 Dilution Buffer  
Immunogen: Electrophorus electricus acetylcholine receptor  
Reported Applications: EAMG induction in rats  
Immunofluorescence  
Formulation: PBS, pH 7.0  
Contains no stabilizers or preservatives  
Endotoxin: <2EU/mg (<0.002EU/μg)  
Determined by LAL gel clotting assay  
Purity: >95%  
Determined by SDS-PAGE  
Sterility: 0.2 μm filtered  
Production: Purified from cell culture supernatant in an animal-free facility  
Purification: Protein G  
RRID: [AB\\_10950153](https://eutils.ncbi.nlm.nih.gov/entrez/eutils/rrid.cgi?db=AB)  
Molecular Weight: 150 kDa

### Description

The mAb35 (TIB-175) antibody reacts with human rat and fish acetylcholine receptor (AChR). The AChR is a heterodimeric receptor consisting of four subunits in a molar stoichiometry  $\alpha 2\beta\gamma\delta$  during the early embryonic stages or after denervation and  $\alpha 2\beta\epsilon\delta$  in the adult form. The mAb35 antibody binds to the major immunogenic region of the  $\alpha$  subunits. It has been shown to passively transfer experimental autoimmune myasthenia gravis.

### Storage

Store at the stock concentration at 4°C . **Do not freeze.**

It is not uncommon for a floccule or precipitate to appear during storage. The floccule is typically buffer salts precipitating out of solution or a small bit of protein aggregation. For information on how to remove floccules or precipitates see our FAQ's at <https://bioxcell.com/faqs>.

### Protocol Information

Since applications vary, each investigator should use the application references as a guide to help estimate the appropriate dose or concentration. The dose or concentration can be further optimized experimentally in a dose response or titration experiment.

### Application References

For a complete list of references, visit [https://bioxcell.com/be0123?bxcs=9k1b3a#tab\\_references](https://bioxcell.com/be0123?bxcs=9k1b3a#tab_references) or scan the QR code below.



---

**Bio X Cell, LLC**  
<https://bioxcell.com>  
+1-866-787-3444  
[customerservice@bioxcell.com](mailto:customerservice@bioxcell.com)

*Conditions: For research use only. Not for use in diagnostic or therapeutic procedures.*

*Not for resale.*

**Bio X Cell, Bio X Cell logo, and all other trademarks are the property of Bio X Cell, LLC © 2025 Bio X Cell, LLC**